




OLEA SPHERE

BREASTSCAPE[®] V1.0

PACKAGE



breastscape[®] v1.0 is the complete breast MR
platform to automate and simplify breast screening.

PRODUCT DESCRIPTION

breastscape® v1.0: intuitive environment for breast MR characterization and follow-up

Automatic subtraction

4D automatic MIP mode

Automatic metrics computation

3D, MIP, MPVR, 3D volume rendering

Lesion follow-up: MR, US, MG DICOM series

Integrated BI-RADS® Atlas report

Multi-screen

PACKAGE EXTENSION

breastscape® v1.0 biopsy extension: MR guided breast interventional procedure planning assistant

Grid template

Identical environment

OTHER AVAILABLE SOLUTIONS

Olea Vision® package: unique solution to easily and productively review MRI and CT images

Body package: complete solution to manage oncology, MSK and metabolic disorders imaging

REQUIREMENTS

System & Software Requirements

Operating System	Windows Desktop (x64): 7 / 8 / 10 (x64) Windows Server (x64): 2008 / 2008 R2 / 2012 / 2012 R2 / 2016 / 2016 R2
Software	Microsoft Visual C++ 2012 Redistributable (x64) Microsoft Visual C++ 2015 Redistributable (x64) PDF Reader

Stand-alone

CPU	Intel® Core / i7 / Xeon processor with 4 cores at 3.0+ GHz or 6 cores at 2.4+ GHz supporting SSE4.2 instructions
Memory	32 GB
GPU	Graphic card with at least the technical specs below: - Memory Type : GDDR5 - Core Clock Speed : 1500MHz - Memory Clock Speed : 8000MHz Compatibles GPU Series: - Nvidia GeForce (Series Kepler, Maxwell, Pascal, Turing) - AMD ATI (Series Radeon Rx)
Storage	150+ GB
Network	1 Gb Ethernet port

Server

CPU	Intel® Core i7 / Xeon processor with 4 cores at 3.0+ GHz or 6 cores at 2.4+ GHz supporting SSE4.2 instructions
Memory	12 GB
Storage	150+ GB
Network	1 Gb Ethernet port

Thick client

CPU	Intel® Core / i7 / Xeon processor with 4 cores at 3.0+ GHz or 6 cores at 2.4+ GHz supporting SSE4.2 instructions
Memory	32 GB
GPU	Graphic card with at least the technical specs below: - Memory Type : GDDR5 - Core Clock Speed : 1500MHz - Memory Clock Speed : 8000MHz Compatibles GPU Series: - Nvidia GeForce (Series Kepler, Maxwell, Pascal, Turing) - AMD ATI (Series Radeon Rx)
Storage	5 GB
Network	1 Gb Ethernet port

Note: Olea Sphere® 3D rendering is not optimal with Nvidia Quadro Series.

OLEA MEDICAL®

ZI Athelia IV - 93, avenue des Sorbiers
13600 La Ciotat - FRANCE
PH +33 4 42 71 24 20 - FX +33 4 42 71 24 27
www.olea-medical.com

Olea Sphere®v3.0 and breastscape® v1.0, medical imaging post-processing software, are medical devices manufactured and marketed by Olea Medical®. These medical devices are reserved for health professionals. These software programs have been designed and manufactured according to the EN ISO 13485 quality management system. Read the instructions in the notice carefully before any use. Instructions for Use are available on <http://www.olea-medical.com/en/>
Manufacturer: Olea Medical® S.A.S. (France). Medical devices Class IIa / Notified body: CE 0459 GMED
BI-RADS ATLAS: a registered trademark of the American College of Radiology (ACR). All rights reserved.

**IMPROVED
DIAGNOSIS
FOR LIFE**